

COMFORT 1100

Functional Emergency Stretcher



- Hydraulic height adjustment
- Trendelenburg and Reverse Trendelenburg
- Backrest adjustment assisted by gas spring
- Viscoelastic cover with antibacterial, air permeable, water-proof polyurethane cover
- X-Ray Translucent, compact laminate mattress platform
- Oxygen tank holder
- Auto regression system
- Lockable and Fold-away PP side rails
- PU pushing bars
- Protective corner bumpers
- Electrostatic powder coated main frame
- Central lockable single castors
- 5th Castor for easy maneuvering
- Height adjustable IV Pole Stand
- ABS Covered sub-frame for easy cleaning and disinfection

Optional accessories:

- Central lockable twin castors (150mm)

TECHNICAL DETAILS

Backrest Angle:	0°-75°
Trendelenburg Angle:	0°-18°
Reverse Trendelenburg Angle:	0°-18°
Castor Diameter:	150mm
Safe Working Load Capacity:	250kg
External Dimensions:	75x210cm
Mattress Dimensions:	60x190cm
Mattress Thickness:	8cm
Height Range: (without mattress)	69-100cm

MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20Şehitkamil Gaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: export@mespa.com.tr

web: www.mespa.com.tr

