

MS 1200



Hydraulic Emergency Stretcher



- Hydraulic height adjustment
- Manual backrest adjustment by gas spring
- Trendelenburg by gas spring
- Reverse Trendelenburg by gas spring
- Easily cleanable, compact laminate mattress platform
- Anti-microbial, water proof, flame retardant, air-permeable mattress
- Collapsible full-length side rails
- Stainless steel undercarriage storage (Optional)
- Lockable, detachable ergonomic pushing handles
- Protective corner bumpers
- IV Pole Sockets
- Height adjustable stainless steel IV pole
- Lockable castors

TECHNICAL DETAILS	
External Dimensions:	72x213cm
Mattress Dimensions:	60x190cm
Mattress Thickness:	8cm
Height (without mattress):	53-81cm
Backrest Angle:	0°-70°
Trendelenburg Angle:	0°-16°
Reverse Trendelenburg Angle:	0°-16°
Castor Diameter:	125mm
Safe Working Load Capacity:	180kg
Weight:	87kg

MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20Şehitkamil Gaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: export@mespa.com.tr

web: www.mespa.com.tr

