

MS 2171



Overbed Table



- Epoxy powder coated metal frame
- Height adjustment by gas spring
- Compact laminate top table
- Lockable castors

TECHNICAL DETAILS	
External Dimensions:	43x80cm
Height Adjustment:	74-98cm
Table Dimensions:	41x80cm
Castor Diameter:	50mm

MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20 ŞehitkamilGaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: export@mespa.com.tr

web: www.mespa.com.tr

