

## Tilt Adjustment Overbed Table



- High pressure laminate (HPL) table with tilting function
- Metal barrier at fixed part of the table
- Gas spring assisted height adjustment
- Thanks to special base, table can be used with beds and cabinets
- Lockable castors

| TECHNICAL DETAILS    |          |
|----------------------|----------|
| External Dimensions: | 39x74cm  |
| Height Adjustment:   | 80-110cm |
| Table Dimensions:    | 39x74cm  |
| Tilting section:     | 39x57cm  |
| Castor Diameter:     | 50mm     |