

VELO BASIC



General Purpose Stretcher



- Manual backrest adjustment
- Height adjustable stainless steel IV pole
- Collapsible, full length side rails with safety lock mechanism
- Anti-microbial, water proof, flame retardant, air-permeable mattress
- Protective corner bumpers
- Detachable ergonomic pushing bars
- IV pole sockets
- Easily cleanable, compact laminate two section mattress platform
- Lockable castors

Optional Accessories:

- Oxygen tube holder for 2liter
- Monitor platform
- Safety belt
- Food table
- Viscoelastic mattress

TECHNICAL DETAILS	
External Dimensions:	72x213cm
Backrest Angle:	0°-75°
Height (without mattress):	60cm
Castor Diameter:	200mm
MattressSurface Dimensions:	60x190cm
Back Section Length:	70cm
Mattress Thickness:	8cm
Safe Working Load Capacity:	250kg

MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20Şehitkamil Gaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: export@mespa.com.tr

web: www.mespa.com.tr

