

MS 1023



Manual ICU and Patient Bed



- Manual adjustment of backrest, footrest, height
- Fowler
- Vascular
- Cardiac
- Collapsible side rails
- 4 section, mattress platform
- Gradually adjustable footrest
- Electrostatic powder coated metal frame
- Lockable and detachable (PP) head and foot boards
- Protective plastic corner bumpers
- IV pole sockets
- Height adjustable stainless steel IV Pole
- Lockable castors




TECHNICAL DETAILS

External Dimensions:	98x214cm
Mattress Platform:	85x195cm
Height(without mattress):	42-75cm
Backrest Angle:	0°-70°
Footrest Angle:	0°-56°
Castor Diameter:	125mm
Safe Working Load Capacity:	250kg
Weight:	105kg

Optional Accessories:

- Lifting pole
- File holder
- Orthopedic traction system
- Trendelenburg and Reverse Trendelenburg

Color Options:



MESPA HEALTHCARE PRODUCTS INC.

2. organize sanayi bolgesi 83216 nolu cadde no:20 Gaziantep / TURKIYE

Tel: +90 444 5 821

Fax: +90 342 337 53 13

e-mail: export@mespa.com.tr

web: www.mespa.com.tr

