



SKE020-1

TABURETE ESCALERAL



SAIKANG / Medical

Healthier
your life

SKE020-1

CARACTERÍSTICAS Y PUNTOS DESTACABLES

Parámetros técnicos

- Size (LxWxH) L490 x W405 x H170-340 mm



Material

Estructura de acero inoxidable con superficie superior de ABS de alta calidad.



Fijo

4 piezas de cubiertas antideslizantes para los pies.



Capa

Escalones dobles.




Embalaje

L520 x W430/170 x H360 mm (1 pieza/caja)
L540 x W430 x H540 mm (2 piezas/caja)





INNOVACIÓN CONSTANTE




Diseño



Investigación



Fabricación



Mantenimiento



Servicio

WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical (Philippines) Company Limited

Manila, The Philippines
TEL: +63 953 044 7864

US Branch

Healthward Technology, INC

17800 CASTLETON ST STE
180.
CITY OF INDUSTRY CA 91748
TEL: +1 (603)416-3063

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 🕒 24/7

✉ export@saikangmedical.com

🌐 www.saikangmedical.com



Company video



General catalogue